

ValProbe (RT) Transport/ Storage Tray

KEEP YOUR DATALOGGERS ORGANIZED AND
SAFE DURING STORAGE AND TRANSPORT



The Logger Tray (V2569) is specifically designed to securely store and transport Kaye ValProbe RT dataloggers, ensuring optimal performance and longevity of your equipment. This tray significantly reduces the risk of damage from accidental drops by providing a secure and organized storage solution. Made from high-quality materials, the tray is robust and

durable, protecting sensitive electronic components from shocks and vibrations while ensuring easy and safe handling during transport.

Utilizing the Logger Tray can not only prevent repair costs and delays but also increase the lifespan and reliability of your dataloggers.


The **ValProbe (RT) Transport/Storage Tray** is an accessory designed to simplify the carrying, storage, and management of the Kaye ValProbe loggers during use. The ValProbe (RT) Transport/Storage Tray can accommodate up to 20 Kaye ValProbe Loggers of any type.

FEATURES

- High-quality, robust design
- Can store up to 20 loggers
- Can be used to securely store the loggers on the shelf
- Ideal for transportation of loggers after programming
- Helpful for sorting non-programmed and pre-

programmed loggers to avoid mistakes if you need to program many loggers

- Made with stainless steel and other materials acceptable for cleanroom environments, including autoclaves
- Suitable for pharma environments, easy to clean with alcohol and other disinfectant sprays
- Rubber feet provide stability when placed on a table or cart
- Rubber feet provides stability when placed on a table or cart
- Size: 300 mm x 210 mm x 95 mm (approximate)
- Includes a unique serial number for easy identification



Kaye Representative contact:

Request a demo:

EUROPE, MIDDLE EAST, AFRICA AND ASIA
 Amphenol Advanced Sensors Germany GmbH
 Sinzheimer Strasse 6
 D-75179 Pforzheim
T: +49 (0) 7231-14 335 0
F: +49 (0) 7231-14335 29
Email: kaye@amphenol-sensors.com

USA/AMERICAS
 Amphenol Thermometrics, Inc.
 967 Windfall Road
 St. Marys, PA 15857
T: +1(814) 834-9140
F: +1(814) 781-7969
Email: kaye-us@amphenol-sensors.com

INDIA
 Amphenol Interconnect India Pvt Ltd.
 Plot no. 6, Survey No.64 | Software Units layout
 MAHAVEER TECHNO PARK
 Hitech City, Madhapur | Hyderabad,
 Telangana – 500081 | **T:** +91 40 33147100
Email: kaye-india@amphenol-sensors.com

CHINA
 Amphenol (Changzhou) Connector Systems
 Co., Ltd, Building 10, Jintong Industrial
 Park, No. 8 Xihu Road, Wujin High-Tech
 Development Zone, Changzhou,
 Jiangsu 213164 | **T:** 0086-519-83055197
Email: kaye-china@amphenol-sensors.com



Warranty and disclaimer: The information mentioned on documents are based on our current tests, knowledge and experience. Because of the effect of possible influences in an application of the product, they do not exempt the user from their own tests, checks and trials. A guarantee of certain properties or a guarantee for the proper suitability of the product for a specific, especially permanent application can not be derived from our data. Liability is therefore excluded to that extent permitted by law. Any proprietary rights of third parties as well as existing laws and regulations must be observed by the recipient of the product on his own responsibility.

© 2024 Amphenol Corporation. All Rights Reserved. Specifications are subject to change without notice. Other company names and product names used in this document are the registered trademarks or trademarks of their respective owners.